

MiBo



Thermoflo IIx Dual

by MIBO MEDICAL GROUP

NEW
BILATERAL
Treatment
Cuts TX time
in HALF





The SOLUTION for patients *suffering* with
CHRONIC DRY EYE

The MiBo Thermoflo is a **therapeutic medical device providing a safe and effective alternative therapy for chronic dry eye.**

- New Dual Hand Piece.
- Quickly convert to Single Hand Piece.
- 5" Touch screen.
- Removable hand pieces.
- Quick startup time.
- Faster temperature recovery.

95% PATIENT
SATISFACTION

Features:

-  **Less Chair Time:** With new Dual Hand Piece cuts chair time in half.
-  **Therapeutic:** Many patients enjoy the soothing massage and spa-like treatment.
-  **No Disposables:** Simply clean the MiBo Thermoflo eye pad between treatments.
-  **Prescribed Therapy Regimen:** With a prescribed therapy regimen the ducts of the Meibomian glands will secrete thinner & clearer lipids for a healthier tear film.

Heater unclips to
convert to a Single
Hand Piece

New Dual
Hand Piece



“If you don't have the MiBo you're not treating dry eye.”

MiBo Medical

Technology for Dry Eye

855-642-6356 | DALLAS, TX
MIBOMEDICALGROUP.COM


All Rights Reserved MIBO Medical Group Rev. 1000A



NOW you can Treat
both eyes at the same
time. Quick turn around
and effective results.

How MIBO THERMOFLO works

Diagram



The MiBo Thermoflo supplies **continuous controlled heat** that is applied to the outer skin of the eye lids along with ultrasound gel for a **gentle massage**.

The heat is absorbed deep into the tissue and **breaks down the hardened lipids** in the meibomian glands.

With a specific **prescribed therapy regimen** the ducts of the meibomian glands will secrete **thinner and clearer lipids** which will allow for a **healthier tear film**.



New C-Elite Eye Pads

Cupped to cover the whole eye. Encapsulates heat into the glands. Requires less movement during treatment.



New Dual Hand Piece

Bilateral treatment reduces time. Adjustable for pupillary distance and orbital comfort.



Benefits for Practice

- Low cost investment
- Affordable treatments
- No disposables
- Tech driven
- Instant relief for patients
- Helpful beyond dry eye!



New User Modes

Now the user has 4 modes to choose from.

Low Mode

-1% Lowest level for comfort and reliability.

Normal Mode

System always defaults to 108.5 degrees F.

High Mode

+1% Higher level for to treat harder lipids.

Ultra Mode

+1.5% Highest level for speciality treatments..

VISIT US ONLINE OR
CALL US TODAY
TO ORDER



Dual Hand Piece
Treatment Modes:
Low 107.5°f, Normal 108.5°f,
High 109.0°f, and Ultra 109.5°f
Weight - 3 lbs
6.25" w x 8.5" d x 4" h
3 Year Limited Warranty
*Single hand piece accessory
sold separately.

ASK The DOCTORS that KNOW about the MiBo Thermoflo

Dr. Harvey Fishman MD PhD Palo Alto, CA

Email: harvfishman@gmail.com

Dr. Thomas Kislak OD Hazleton, PA

Email: foreeyes@ptd.net

Dr. James Lewis MD Philadelphia, PA

Email: jslewis@jameslewismd.com

* These doctors are not paid employees of MIBO Medical Group.

MIBO Medical Group respects these doctors for taking their time to tell others about our MiBoThermoFlo and we hope you do the same.

MIBO Medical Group is a privately owned company that started in October 2012. Our products are made in Dallas, TX with strict ISO 13485 manufacturing process. The MiBo is a FDA class II with exempt therapeutic medical device.

MiBo Medical

Technology for Dry Eye

855-642-6356 | DALLAS, TX
MIBOMEDICALGROUP.COM

All Rights Reserved MIBO Medical Group Rev. 1000A